











BRUDER 5 MM Ø LAPAROSCOPIC INSTRUMENTS

**ADDITIONAL HANDLE STYLES
SPECIFY NUMBER WHEN ORDERING**

Handles available for all Premier 5 mm Rigid, Rotatable and Modular style instruments

3B.) ADDITIONAL HANDLE STYLES 5-12:		**ADDITIONAL CHARGE FOR HANDLES 5-12**	
			
<p>TYPE 5: LIGHTWEIGHT CARBON FIBER WITH RF-POST</p> <p>FINGER TANG FOR EXTRA CONTROL, ERGONOMICAL MOLDED DESIGN, FULLY AUTOCLAVABLE</p>	<p>TYPE 6: LIGHTWEIGHT CARBON FIBER WITH RF-POST, PRECISE ENGAGEMENT RATCHET</p> <p>FINGER TANG FOR EXTRA CONTROL, ERGONOMICAL MOLDED DESIGN, FULLY AUTOCLAVABLE</p>	<p>TYPE 7: LIGHTWEIGHT CARBON FIBER WITH RF-POST, ON-DEMAND ON/OFF RATCHET</p> <p>FINGER TANG FOR EXTRA CONTROL, ERGONOMICAL MOLDED DESIGN, FULLY AUTOCLAVABLE</p>	<p>TYPE 8: LIGHTWEIGHT CARBON FIBER WITH RF-POST, FRENCH ON/OFF RATCHET</p> <p>ERGONOMICAL MOLDED DESIGN, FULLY AUTOCLAVABLE</p>
			
<p>TYPE 9: LIGHTWEIGHT CARBON FIBER WITH RF-POST</p> <p>EXTRA LARGE FINGER LOOPS FOR LARGER HANDS AND EXTRA POWER, ERGONOMICAL MOLDED DESIGN, FULLY AUTOCLAVABLE</p>	<p>TYPE 10: LIGHTWEIGHT CARBON FIBER WITH RF-POST, PRECISE ENGAGEMENT RATCHET</p> <p>EXTRA LARGE FINGER LOOPS FOR LARGER HANDS AND EXTRA POWER, ERGONOMICAL MOLDED DESIGN, FULLY AUTOCLAVABLE</p>	<p>TYPE 11: STAINLESS STEEL AXIAL HANDLE</p> <p>ERGONOMICAL DESIGN FOR EXTRA CONTROL, BALANCE & POWER W/ FLUSH PORT <i>(Not available on modular style instrument. Additional cost for rotatable style instrument)</i></p>	<p>TYPE 12: STAINLESS STEEL AXIAL HANDLE W/ ON/OFF RATCHET LOCK</p> <p>ERGONOMICAL DESIGN FOR EXTRA CONTROL, BALANCE & POWER W/ FLUSH PORT <i>(Not available on modular style instrument. Additional cost for rotatable style instrument)</i></p>

NEW FROM BRUDER...

Always a sharp, precise cut with the BRUDER **COST-SAVING METZENBAUM Laparoscopic Scissors** Disposable sterile tip and reusable insulated stainless steel handle



Sterile Tip
(Box of 10 individually packaged)

ORDERING INFORMATION:

List No. ARO60001	CURVED METZENBAUM TIP-STERILE (PKG. OF 10)
List No. ARO60002	STRAIGHT METZENBAUM TIP-STERILE (PKG. OF 10)
List No. ARO60510	5 MM DIA., 33 CM INSULATED STAINLESS STEEL HANDLE AND SHAFT, ROTATABLE 1-PIECE WITH FLUSH PORT
List No. ARO61510	5 MM DIA., 45 CM INSULATED STAINLESS STEEL HANDLE AND SHAFT, ROTATABLE 1-PIECE WITH FLUSH PORT