



# THE ORIGINAL KLEPPINGER BIPOLAR LAPAROSCOPIC FORCEPS

## AVAILABLE WITH REUSABLE OPTIONAL TIP INNER TONGS

**REUSABLE REPLACEMENT TONGS FOR KLEPPINGER TYPE FORCEPS:** (All of our inner tongs are interchangeable with the original Wolf style.)



**FLAT TIP (F)**  
CONVENTIONAL KLEPPINGER  
STYLE



**STANDARD TIP (S)**  
FENESTRATED W/ WAVE  
FOR SECURE, SLIP-FREE GRIP



**FENESTRATED TIP (T)**  
W/ TOOTH  
FOR SECURE, SLIP-FREE GRIP



**MICRO TIP (M)**  
FOR PRECISE,  
PIN-POINT COAGULATION

LIST NO. <b>ARO205101</b>	32 CM - FLAT TIP (ORIGINAL KLEPPINGER) INNER TONG - F
LIST NO. <b>ARO105101</b>	45 CM - FLAT TIP (ORIGINAL KLEPPINGER) INNER TONG - F
LIST NO. <b>ARO205102</b>	32 CM - STANDARD TIP (FENESTRATED WITH WIRE) - INNER TONG - S
LIST NO. <b>ARO105102</b>	45 CM - STANDARD TIP (FENESTRATED WITH WIRE) - INNER TONG - S
LIST NO. <b>ARO205103</b>	32 CM - MICRO TIP (NARROW SHORT TIP) - INNER TONG - M
LIST NO. <b>ARO105103</b>	45 CM - MICRO TIP (NARROW SHORT TIP) - INNER TONG - M
LIST NO. <b>ARO205104</b>	32 CM - FENESTRATED TIP (FENESTRATED WITH TOOTH) - INNER TONG - T
LIST NO. <b>ARO105104</b>	45 CM - FENESTRATED TIP (FENESTRATED WITH TOOTH) - INNER TONG - T

**Note:** Both the syringe and scissor handle Kleppinger forceps are available with the optional inner tongs. Use the first number from the front page for the instrument and add the letter of your choice for the optional insert S, T, M.

**BRUDER** also has reusable replacement tongs available for Storz, ACMI, Piling etc and other brands of bipolar laparoscopic forceps.

### BIPOLAR CABLES:

BRUDER RF-Cables are fully autoclavable, 10 feet long



BRUDER/WOLF  
KLEPPINGER PLUG &  
BRUDER/WOLF  
GENERATOR PLUG

LIST NO. **ARO00523**



BRUDER/WOLF  
KLEPPINGER PLUG &  
DUAL-BANANA  
GENERATOR PLUG

LIST NO. **ARO00524**

### CLEANING BRUSHES:

BRUDER cleaning brushes are reusable



LIST NO. <b>ARO00329</b>	CLEANING BRUSHES FOR INNER SHEATH, 19" LONG, (PER DOZEN)
LIST NO. <b>ARO00331</b>	CLEANING BRUSHES FOR OUTER SHEATH, 19" LONG, (PER DOZEN)
LIST NO. <b>ARO00343</b>	FLAT BRUSH WITH HANDLE FOR OUTSIDE SURFACE

### LBC 50-P DIGITAL

LIST NO. **ARO5671**  
SEE BROCHURE FOR  
INFORMATION

#### FEATURES:

- Audible and Visual coagulation monitor, "End Tone" and automatic output power
- Digital wattage display for easy-to-reproduce settings



### ARFM RF-METER

LIST NO. **ARO5675** OR  
**ARO5680** FOR USE WITH  
OTHER GENERATORS



For safe and efficient tubular sterilization **BRUDER** offers the **LBC 50-P Bipolar Generator**, with both visual display and audio signal for coagulation monitoring, along with the **ARFM Bipolar Laparoscopic Coagulation Monitor** for use in conjunction with bipolar generators.

For complete Laparoscopic Instrument information refer to our applicable data sheets

Save and buy direct from the manufacturer.  
We design, manufacture & sell the tools the surgeons use.

**BRUDER** specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



**ARROW MEDICAL SUPPLY INC**  
Telephone Number: (847) 367-7771  
Web Site: [www.arrowmedical.com](http://www.arrowmedical.com)

28085 N Ashley Circle Bldg 105

Libertyville IL 600486  
Fax Number: (847)367-7751  
Email: [info@arrowmedical.com](mailto:info@arrowmedical.com)