




## BRUDER LAPAROSCOPY

### BRUDER RETRACTORS FOR MINIMAL INVASIVE SURGERY


**New from  
BRUDER...**

#### ARTICULATING FAN RETRACTOR

- All stainless steel construction with flush port
- 60° articulating folding fan blades



LIST NO. <b>ARO50515</b>	<b>5 MM Ø, 34 CM ARTICULATING FAN RETRACTOR</b> RETRACTOR / ELEVATOR WITH 3 FLAT, 60 DEGREE ARTICULATING AND FOLDING BLADES, (60 MM MAX WIDTH), ADJUSTABLE HANDLE WITH FLUSHPORT, STAINLESS STEEL
LIST NO. <b>ARO51015</b>	<b>10 MM Ø, 34 CM ARTICULATING FAN RETRACTOR</b> RETRACTOR / ELEVATOR WITH 5 FLAT, 60 DEGREE ARTICULATING AND FOLDING BLADES, (90 MM MAX WIDTH), ADJUSTABLE HANDLE WITH FLUSHPORT, STAINLESS STEEL
LIST NO. <b>ARO52515</b>	<b>5 MM Ø, 45 CM ARTICULATING FAN RETRACTOR</b> RETRACTOR / ELEVATOR WITH 3 FLAT, 60 DEGREE ARTICULATING AND FOLDING BLADES, (60 MM MAX WIDTH), ADJUSTABLE HANDLE WITH FLUSHPORT, STAINLESS STEEL
LIST NO. <b>ARO51515</b>	<b>10 MM Ø, 45 CM ARTICULATING FAN RETRACTOR</b> RETRACTOR / ELEVATOR WITH 5 FLAT, 60 DEGREE ARTICULATING AND FOLDING BLADES, (90 MM MAX WIDTH), ADJUSTABLE HANDLE WITH FLUSHPORT, STAINLESS STEEL



LIST NO. <b>ARO50550</b>	<b>5 MM Ø, 33 CM RETRACTOR/ELEVATOR</b> CURVED, ROTATABLE, TRIANGULAR KITE-SHAPED BLADE, LOCKING HANDLE
LIST NO. <b>ARO51550</b>	<b>5 MM Ø, 45 CM RETRACTOR/ELEVATOR</b> CURVED, ROTATABLE, TRIANGULAR KITE-SHAPED BLADE, LOCKING HANDLE
LIST NO. <b>ARO50551</b>	<b>5 MM Ø, 33 CM RETRACTOR/ELEVATOR</b> STRAIGHT, ROTATABLE, TRIANGULAR KITE-SHAPED BLADE, LOCKING HANDLE
LIST NO. <b>ARO51551</b>	<b>5 MM Ø, 45 CM RETRACTOR/ELEVATOR</b> STRAIGHT, ROTATABLE, TRIANGULAR KITE-SHAPED BLADE, LOCKING HANDLE
LIST NO. <b>ARO51050</b>	<b>10 MM Ø, 33 CM RETRACTOR/ELEVATOR</b> STRAIGHT, ROTATABLE, TRIANGULAR KITE-SHAPED BLADE, LOCKING HANDLE
LIST NO. <b>ARO51150</b>	<b>10 MM Ø, 45 CM RETRACTOR/ELEVATOR</b> STRAIGHT, ROTATABLE, TRIANGULAR KITE-SHAPED BLADE, LOCKING HANDLE

## BRUDER RETRACTORS FOR MINIMAL INVASIVE SURGERY



<b>LIST NO.</b> <b>ARO50539</b>	<b>5 MM Ø, 33 CM LEVER/ELEVATOR</b> RETRACTABLE, STRAIGHT, 4-FINGER BLADES, ROTATABLE INSERT
<b>LIST NO.</b> <b>ARO51539</b>	<b>5 MM Ø, 45 CM LEVER/ELEVATOR</b> RETRACTABLE, STRAIGHT, 4-FINGER BLADES, ROTATABLE INSERT




<b>LIST NO.</b> <b>ARO20551</b>	<b>5 MM Ø, 33 CM ELEVATOR/PROBE</b> STRAIGHT, TEFLON TIP, ACCORDING TO DR. VALLE
------------------------------------	---



<b>LIST NO.</b> <b>ARO50540</b>	<b>5 MM Ø, 33 CM RETRACTOR/ELEVATOR</b> 90 DEGREE ARTICULATING TEFLON TIP, ACCORDING TO DR. VALLE, WITH FLUSH PORT
------------------------------------	---



<b>LIST NO.</b> <b>ARO50548</b>	<b>5 MM Ø, 33 CM RETRACTOR</b> 90 DEGREE ARTICULATING FLAT BLADE, WITH FLUSH PORT
<b>LIST NO.</b> <b>ARO51548</b>	<b>5 MM Ø, 45 CM RETRACTOR</b> 90 DEGREE ARTICULATING FLAT BLADE, WITH FLUSH PORT
<b>LIST NO.</b> <b>ARO51087</b>	<b>10 MM Ø, 33 CM RETRACTOR</b> 90 DEGREE ARTICULATING FLAT BLADE, WITH FLUSH PORT
<b>LIST NO.</b> <b>ARO51187</b>	<b>10 MM Ø, 45 CM RETRACTOR</b> 90 DEGREE ARTICULATING FLAT BLADE, WITH FLUSH PORT
<b>LIST NO.</b> <b>ARO51086</b>	<b>10 MM Ø, 33 CM RETRACTOR</b> 90 DEGREE ARTICULATING "CUSHING"-TYPE BLADE, WITH FLUSH PORT
<b>LIST NO.</b> <b>ARO51186</b>	<b>10 MM Ø, 45 CM RETRACTOR</b> 90 DEGREE ARTICULATING "CUSHING"-TYPE BLADE, WITH FLUSH PORT



Save and buy direct from the manufacturer.  
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



**ARROW MEDICAL SUPPLY INC**  
Telephone Number: (847) 367-7771  
Web Site: [www.arrowmedical.com](http://www.arrowmedical.com)

28085 N Ashley Circle Bldg 105

Libertyville IL 600486  
Fax Number: (847)367-7751  
Email: [info@arrowmedical.com](mailto:info@arrowmedical.com)