











## BRUDER LAPAROSCOPY

### 5 MM Ø, LAPAROSCOPIC MONOPOLAR ELECTRODES FOR PRECISE DISSECTING AND COAGULATING

Our line of monopolar electrodes are completely reusable and are manufactured to the highest standards. We offer a wide variety of tip designs, both solid shaft and suction models with unique designed handles for ergonomic ease of handling and electrostatic powder coating insulation for maximum safety and micro polished tip for minimum eschar buildup.

### LAPAROSCOPIC MONOPOLAR ELECTRODES













Maximum Insulated Tip Design	Shaft Length		
		Side Mount RF Connector	Rear Mounted RF Connector
 "J" Hook	33 cm	List No. <b>ARO20538</b>	List No. <b>ARO20581</b>
	45 cm	List No. <b>ARO10533</b>	List No. <b>ARO10581</b>
	CUSTOM	CALL ARROW	CALL ARROW
 "L" Hook	33 cm	List No. <b>ARO20578</b>	List No. <b>ARO20582</b>
	45 cm	List No. <b>ARO10535</b>	List No. <b>ARO10582</b>
	CUSTOM	CALL ARROW	CALL ARROW
 Spatula	33 cm	List No. <b>ARO20556</b>	List No. <b>ARO20583</b>
	45 cm	List No. <b>ARO10546</b>	List No. <b>ARO10583</b>
	CUSTOM	CALL ARROW	CALL ARROW
 Cone	33 cm	List No. <b>ARO20536</b>	List No. <b>ARO20588</b>
	45 cm	List No. <b>ARO10531</b>	List No. <b>ARO10588</b>
	CUSTOM	CALL ARROW	CALL ARROW
 Ball	33 cm	List No. <b>ARO20537</b>	List No. <b>ARO20589</b>
	45 cm	List No. <b>ARO10532</b>	List No. <b>ARO10589</b>
	CUSTOM	CALL ARROW	CALL ARROW
 Knife	33 cm	List No. <b>ARO20535</b>	List No. <b>ARO20590</b>
	45 cm	List No. <b>ARO10530</b>	List No. <b>ARO10590</b>
	CUSTOM	CALL ARROW	CALL ARROW

**Needle electrodes are featured on our Needle Electrode Data Sheet.  
All instruments are available in 3 mm Ø - See Applicable Data Sheet.**

# LAPAROSCOPIC MONOPOLAR ELECTRODES WITH SUCTION PORT



Maximum Insulated Tip Design	Shaft Length h			
		Luer Lock Connection & Side RF Post	Stopcock Connection & Side RF Post	Trumpet Valve Connection & Side RF Post
 "J" Hook	33 cm	List No. <b>ARO20557L</b>	List No. <b>ARO20557</b>	List No. <b>ARO20563</b>
	45 cm	List No. <b>ARO10548L</b>	List No. <b>ARO10548</b>	List No. <b>ARO10554</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW
 "L" Hook	33 cm	List No. <b>ARO20572L</b>	List No. <b>ARO20572</b>	List No. <b>ARO20573</b>
	45 cm	List No. <b>ARO10556L</b>	List No. <b>ARO10556</b>	List No. <b>ARO10557</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW
 Spatula	33 cm	List No. <b>ARO20558L</b>	List No. <b>ARO20558</b>	List No. <b>ARO20564</b>
	45 cm	List No. <b>ARO10553L</b>	List No. <b>ARO10553</b>	List No. <b>ARO10550</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW
 Cone	33 cm	List No. <b>ARO20577L</b>	List No. <b>ARO20577</b>	List No. <b>ARO20584</b>
	45 cm	List No. <b>ARO10558L</b>	List No. <b>ARO10558</b>	List No. <b>ARO10560</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW
 Sling	33 cm	List No. <b>ARO20566L</b>	List No. <b>ARO20566</b>	List No. <b>ARO20568</b>
	45 cm	List No. <b>ARO10549L</b>	List No. <b>ARO10549</b>	List No. <b>ARO10552</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW
 Knife	33 cm	List No. <b>ARO20565L</b>	List No. <b>ARO20565</b>	List No. <b>ARO20567</b>
	45 cm	List No. <b>ARO10547L</b>	List No. <b>ARO10547</b>	List No. <b>ARO10551</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW
 Suction Tip	33 cm	List No. <b>ARO20528L</b>	List No. <b>ARO20527</b>	List No. <b>ARO20526</b>
	45 cm	List No. <b>ARO10528L</b>	List No. <b>ARO10528</b>	List No. <b>ARO10528T</b>
	CUSTOM	CALL ARROW	CALL ARROW	CALL ARROW

## REUSABLE MONOPOLAR RF-CABLES FOR USE WITH ELECTRODES:

List No. <b>ARO00531</b>	Monopolar RF Cable with Female Banana (4 mm Ø) Instrument Connector and Male Banana (4 mm Ø) unit connector
List No. <b>ARO930</b>	Monopolar RF Cable with Female Banana (4 mm Ø) Instrument Connector and Male Bovie (8 mm Ø) unit connector

Save and buy direct from the manufacturer.  
We design, manufacture & sell the tools the surgeons use.

**BRUDER** specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



**ARROW MEDICAL SUPPLY INC**  
Telephone Number: (847) 367-7771  
Web Site: [www.arrowmedical.com](http://www.arrowmedical.com)

28085 N Ashley Circle Bldg 105

Libertyville IL 600486  
Fax Number: (847)367-7751  
Email: [info@arrowmedical.com](mailto:info@arrowmedical.com)