













BRUDER LAPAROSCOPY

10 MM Ø DISSECTING, GRASPING, BIOPSY FORCEPS AND SCISSORS











BRUDER, one of the leaders in design and manufacturing of laparoscopic instruments offer the most complete line of 10 mm Ø instruments. Our new line of functional style A, B and C instruments provide the ultimate in versatility, they are completely reusable and competitively priced.

BRUDER Multifunction Style Laparoscopic Instruments	
A	Rigid, Non-Take-Apart Instrument with flush port – Choose Handle
B	Rotatable, Non-Take-Apart with flush port – Choose Handle
C	Rotatable, Modular – Take-A-Part, 3-piece Instrument – Choose Handle





When ordering, add the letter of the functional style of instrument and handle type number to the list number below.

	Strong Grasper/Dissector , insulated shaft, suggested insulated handle with RF-connector (Recommend Type A Handle)
	33 cm List No. ARO91002 45 cm List No. ARO910102
	Maryland Grasper/Dissector , insulated shaft, suggested insulated handle with RF-connector (Recommend Type A Handle)
	33 cm List No. ARO91006 45 cm List No. ARO910106
	Mixer Grasper/Dissector , right angle, insulated shaft, suggested insulated handle with RF-connector (Recommend Type A Handle)
	33 cm List No. ARO91007 45 cm List No. ARO910107
	Mixer Grasper/Dissector , 45° angle, insulated shaft, suggested insulated handle with RF-connector (Recommend Type A Handle)
	33 cm List No. ARO91008 45 cm List No. ARO910108
	Babcock Grasper , fenestrated tip, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)
	33 cm List No. ARO91016 45 cm List No. ARO910116
	US Babcock Grasper , spring tension, atraumatic, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)
	33 cm List No. ARO91017 45 cm List No. ARO910117
	Allis Grasper , fine axial serrations, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)
	33 cm List No. ARO91021 45 cm List No. ARO910121
	DeBakey Grasper , 40 mm jaw, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)
	33 cm List No. ARO91022 45 cm List No. ARO910122
	DeBakey Grasper , 60 mm jaw, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)
	33 cm List No. ARO91023 45 cm List No. ARO910123
	Kocher Grasper , 40 mm jaw, axial serrations, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)
	33 cm List No. ARO91024 45 cm List No. ARO910124

10 MM Ø DISSECTING, GRASPING, BIOPSY FORCEPS AND SCISSORS

	Kocher Grasper , 60 mm jaw, axial serrations, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)			
	33 cm	List No. ARO91025	45 cm	List No. ARO910125
	Duval Grasper , 40 mm jaw, axial serrations, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)			
	33 cm	List No. ARO91026	45 cm	List No. ARO910126
	Claw Forceps , 2 x 3 teeth, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)			
	33 cm	List No. ARO91036	45 cm	List No. ARO910136
	Spoon Forceps , standard tip, insulated shaft, suggested insulated ratchet handle with RF Post (Recommend Type B Handle)			
	33 cm	List No. ARO91037	45 cm	List No. ARO910137
	Spoon Forceps , long tip, insulated shaft, suggested insulated ratchet handle with RF Post (Recommend Type B Handle)			
	33 cm	List No. ARO91038	45 cm	List No. ARO910138
	Schroeder Tenaculum Forceps , single tooth, insulated shaft, suggested non-insulated ratchet handle (Recommend Type B Handle)			
	33 cm	List No. ARO91045	45 cm	List No. ARO910145
	Strong Atraumatic Clamp , 40 mm straight jaw, insulated shaft, suggested insulated ratchet handle with RF Post (Recommend Type B Handle)			
	33 cm	List No. ARO91046	45 cm	List No. ARO910146
	Strong Atraumatic Clamp , 40 mm curved jaw, insulated shaft, suggested insulated ratchet handle with RF post (Recommend Type B Handle)			
	33 cm	List No. ARO91047	45 cm	List No. ARO910147
	Straight Scissors , single action, insulated shaft, suggested handle with RF-connector (Recommend Type A Handle)			
	33 cm	List No. ARO91064	45 cm	List No. ARO910164
	Curved Scissors , double action, insulated shaft, suggested insulated handle with RF-connector (Recommend Type A Handle)			
	33 cm	List No. ARO91066	45 cm	List No. ARO910166

BRUDER Handles are made of stainless steel, insulated as required by electrostatic powder coating
Standard Handles shown in the text can be substituted for the handles of your choice at no extra charge.

			
Type A Insulated, RF- Connector	Type 2 Insulated, RF-Connector, Ratchet Handle	Type 3 Un-insulated	Type B Un-insulated, Ratchet Handle

Save and buy direct from the manufacturer.
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



ARROW MEDICAL SUPPLY INC
Telephone Number: (847) 367-7771
Web Site: www.arrowmedical.com

28085 N Ashley Circle Bldg 105

Libertyville IL 600486
Fax Number: (847)367-7751
Email: info@arrowmedical.com