

**Serving the Medical
Profession Since
1990**



847-367-7771

**BRUDER SURGICAL
INSTRUMENTATION**


*By
ARROW MEDICAL SUPPLY INC*

BRUDER LAPAROSCOPY

TULANDI KNOT PUSHER

Designed in conjunction with Togas Tulandi, M.D.

The **TULANDI KNOT PUSHER** has a specially designed tip for performing extracorporeal knotting with speed, ease of execution and reliability.



LIST NO. ARO20561	TULANDI KNOT PUSHER STAINLESS STEEL 5 MM Ø, 33 CM LONG
--------------------------	--

In addition the suture manipulator (knot pusher) can be used to form a loop and the rounded paddle-shaped distal end makes it possible to use the instrument as a normal palpation probe.

Like all BRUDER Laparoscopic instruments, this instrument is handcrafted by expert technicians utilizing the finest grade Surgical Stainless Steel; completely reusable & competitively priced.

Save and buy direct from the manufacturer.
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



ARROW MEDICAL SUPPLY INC
Telephone Number: (847) 367-7771
Web Site: www.arrowmedical.com

28085 N Ashley Circle Bldg 105

Libertyville IL 600486
Fax Number: (847)367-7751
Email: info@arrowmedical.com