



BRUDER LAPAROSCOPY

OPERATING LAPAROSCOPES WITH SUPERIOR OPTICS

BRUDER has been a pacesetter in the design and manufacturing of fine quality laparoscopic equipment and instrumentation for over 20 years. **BRUDER** continues to lead the way with quality built laparoscopes to satisfy the most demanding criteria.

All **BRUDER** laparoscopes are hand crafted by journeyman optical instrument makers utilizing the highest grade materials. They are completely manufactured in our facility in Addison, IL, USA.

BRUDER AUTOCLAVABLE ROD LENS LAPAROSCOPES

Our autoclavable rod lens laparoscopes represent the very latest in technology in the design and fabrication of rigid endoscopes. The expanded field of view and light transmission characteristics are unsurpassed by any other rod lens optic on the market. Our new rod lens laparoscopes have their entire critical seams laser welded to exact tolerances. The objective and eyepiece lenses are secured with a unique adhesive and coating material. Such innovations make these laparoscopes much more durable and fully autoclavable.



List No. ARO20501



List No. ARO1030B



List No. ARO51501

LIST NO. ARO20500	LAPAROSCOPE: 0 DEGREE VIEWING, 5 MM DIAMETER, 260 MM LONG, AUTOCLAVABLE
LIST NO. ARO20501	LAPAROSCOPE: 30 DEGREE VIEWING, 5 MM DIAMETER, 260 MM LONG, AUTOCLAVABLE
LIST NO. ARO5000B	LAPAROSCOPE: 0 DEGREE VIEWING, 5 MM DIAMETER, 300 MM LONG, AUTOCLAVABLE
LIST NO. ARO5030B30	LAPAROSCOPE: 30 DEGREE VIEWING, 5 MM DIAMETER, 300 MM LONG, AUTOCLAVABLE
LIST NO. ARO1000B	LAPAROSCOPE: 0 DEGREE VIEWING, 10 MM DIAMETER, 300 MM LONG, AUTOCLAVABLE
LIST NO. ARO1030B	LAPAROSCOPE: 30 DEGREE VIEWING, 10 MM DIAMETER, 300 MM LONG, AUTOCLAVABLE
LIST NO. ARO1030B45	LAPAROSCOPE: 45 DEGREE VIEWING, 10 MM DIAMETER, 300 MM LONG, AUTOCLAVABLE
LIST NO. ARO51500	LAPAROSCOPE: 0 DEGREE VIEWING, 10 MM DIAMETER, 400 MM LONG (FOR BARIATRIC SURGERY), AUTOCLAVABLE
LIST NO. ARO51501	LAPAROSCOPE: 30 DEGREE VIEWING, 10 MM DIAMETER, 400 MM LONG (FOR BARIATRIC SURGERY), AUTOCLAVABLE

BRUDER OPERATING LAPAROSCOPES WITH SUPERIOR OPTICS

- Fully Autoclavable
- German Optic Rod Lens System
- Hand crafted by journeyman optical instrument makers utilizing the highest grade stainless steel and polymers
- ACMI, Wolf, Storz / Olympus Light Cable connectors
- We offer the widest range of operating instrument channels
- Made with pride in the USA by BRUDER



List No. ARO11003

LIST NO. ARO11001	BRUDER OPERATING LAPAROSCOPE 0 DEGREE, 10 MM Ø, RIGHT ANGLE EYE PIECE, 3 MM INSTRUMENT CHANNEL
LIST NO. ARO11002	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 10 MM Ø, 45 DEGREE ANGLE EYE PIECE, 3 MM INSTRUMENT CHANNEL
LIST NO. ARO11003	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 10 MM Ø, RIGHT ANGLE EYE PIECE, 5.5 MM INSTRUMENT CHANNEL
LIST NO. ARO11004	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 11 MM Ø, 45 DEGREE ANGLE EYE PIECE, 6 MM INSTRUMENT CHANNEL
LIST NO. ARO11005	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 11 MM Ø, RIGHT ANGLE EYE PIECE, 6 MM INSTRUMENT CHANNEL
LIST NO. ARO11006	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 12 MM Ø, RIGHT ANGLE EYE PIECE, 7 MM INSTRUMENT CHANNEL
LIST NO. ARO11007	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 12 MM Ø, 45 DEGREE ANGLE EYE PIECE, 7 MM INSTRUMENT CHANNEL
LIST NO. ARO11008	BRUDER OPERATING LAPAROSCOPE 0 DEGREE VIEWING, 10 MM Ø, 45 DEGREE ANGLE EYE PIECE, 5.5 MM INSTRUMENT CHANNEL

Save and buy direct from the manufacturer.
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



ARROW MEDICAL SUPPLY INC
Telephone Number: (847) 367-7771
Web Site: www.arrowmedical.com

28085 N Ashley Circle Bldg 105

Libertyville IL 600486
Fax Number: (847)367-7751
Email: info@arrowmedical.com