

Serving the Medical
Profession Since
1990



**BRUDER SURGICAL
INSTRUMENTATION**

By
ARROW MEDICAL SUPPLY INC




BRUDER BARIATRIC LAPAROSCOPY

BARIATRIC SURGERY WITH THE RETRACT-ROBOT

(MADE WITH PRIDE IN THE USA)

THE NEVER TIRING, TOTALLY OBEDIENT SURGICAL ASSISTANT

The most versatile tool for the Operating Room and the most quality cost efficient operating room purchase you can make.



*Bariatric surgery performed by Dr. Mark Colquitt, M.D. of Blount Memorial Medical Center, Maryville, TN, utilizing the **Retract-Robot** (list # ARO08801, with 26'' long shaft) and the Nathanson's style liver retractor (list #ARO42000 or #ARO43000).*

The Retractable-Robot offers greater freedom and ease of use for the surgeon.

- The Retractable-Robot is attached over the surgical drape of the operating table with our sterile clamp (list no. ARO09056) to the side rail for greater room, operative convenience and sterility
- Retractor placement and pressure can be controlled for optimal exposure
- Stable and constant traction, as required by the surgeon
- Less clutter in the operating field and greater exposure equals improved performance
- Minimizes staff requirements


BARIATRIC SURGERY WITH THE RETRACT-ROBOT

BRUDER presents...

The ideal instrument for Bariatric Surgery for liver retraction/elevation using the RETRACT-ROBOT and Nathanson Style Liver Retractor or Kite Retractor.



Retract-Robot holding a Nathanson Liver Retractor & Goose Neck Retractor holding Laparoscope



Retract-Robot holding a Nathanson Liver Retractor.



Retract-Robot holding an Articulating Fan Retractor.

The RETRACT-ROBOT:

The never-tiring, totally obedient surgical assistant for bariatric surgery & retractor used in conjunction

ORDERING INFORMATION:

| | |
|----------|---|
| ARO08801 | RETRACT-ROBOT WITH EXTRA LONG ARM |
| ARO09053 | STERIL CLAMP: AUTOCLAVABLE CLAMP FOR MOUNTING TO THE TABLE ABOVE THE SURGICAL DRAPE |
| ARO42000 | NATHANSON STYLE LIVER RETRACTORS, 6.5MM DIA., SET OF 3 (SMALL, MEDIUM, LARGE) |
| ARO43000 | NATHANSON STYLE LIVER RETRACTORS, 5MM DIA., SET OF 3 (SMALL, MEDIUM, LARGE) |
| ARO51515 | ARTICULATING FAN RETRACTOR/ELEVATOR, 5MM DIA., 37CM |
| ARO51015 | ARTICULATING FAN RETRACTOR/ELEVATOR, 10MM DIA., 37CM |
| ARO51150 | CURVED KITE-SHAPED RETRACTOR, 5MM DIA., 45CM |
| ARO50550 | STRAIGHT KITE-SHAPED RETRACTOR, 10MM DIA., 45CM |

Save and buy direct from the manufacturer.
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



ARROW MEDICAL SUPPLY INC
Telephone Number: (847) 367-7771
Web Site: www.arrowmedical.com

28085 N Ashley Circle Bldg 105

Libertyville IL 600486

Fax Number: (847)367-7751

Email: info@arrowmedical.com