





BRUDER LAPAROSCOPY

BRUDER MONOPOLAR NEEDLE ELECTRODES FOR LAPAROSCOPIC SURGERY

EXCLUSIVE BRUDER RETRACTABLE NEEDLE ELECTRODES FOR PRECISE LAPAROSCOPIC CUTTING & COAGULATION



LIST NO. ARO20553	30 CM INSULATED CUTTING ELECTRODE, RETRACTABLE NEEDLE, ACCORDING TO DR. VALLE, 5 MM DIAMETER SHAFT WITH 3 MM TIP FOR BETTER VISION, SET SCREW ADJUSTMENT (FIGURE 1)
LIST NO. ARO20554	REPLACEMENT STYLET, FOR #ARO20553
LIST NO. ARO20569	30 CM INSULATED CUTTING ELECTRODE, RETRACTABLE NEEDLE, WITH IRRIGATION BY TRUMPET VALVE, ACCORDING TO DR. VALLE, 5 MM DIAMETER SHAFT, SET SCREW ADJUSTMENT
LIST NO. ARO20570	IRRIGATION SLEEVE ONLY, FOR #ARO20569
LIST NO. ARO205110	30 CM INSULATED CUTTING ELECTRODE, RETRACTABLE NEEDLE 5 MM DIAMETER SHAFT WITH 3 MM TIP FOR BETTER VISION, SPRING LOADED PUSH BUTTON (FIGURE 2)



LIST NO. ARO20548	30 CM INSULATED CUTTING ELECTRODE, WITH CHUCK FOR MICRO NEEDLES, ACCORDING TO DR. VALLE, 5 MM DIAMETER SHAFT (FIGURE 3)
LIST NO. ARO00533	REPLACEMENT 150 MICRON NEEDLE FOR #ARO20548
LIST NO. ARO00534	REPLACEMENT 500 MICRON NEEDLE FOR #ARO20548

**THE SAME MONOPOLAR NEEDLE ELECTRODES ARE AVAILABLE IN 44 CM.
PLEASE REQUEST ORDER NUMBER AND PRICE.**

Save and buy direct from the manufacturer.
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



ARROW MEDICAL SUPPLY INC
Telephone Number: (847) 367-7771
Web Site: www.arrowmedical.com

28085 N Ashley Circle Bldg 105

Libertyville IL 600486
Fax Number: (847)367-7751
Email: info@arrowmedical.com