



BRUDER LAPAROSCOPY

LAPAROSCOPIC INSTRUMENTATION FOR BARIATRIC SURGERY



BRUDER, one of the leaders in Endoscopic instruments for over two decades, in conjunction with leading surgeons, has developed an extensive line of bariatric laparoscopic instruments to meet the growing demand for surgical treatment of morbid obesity.

With the recent advances in laparoscopy, the laparoscopic approach when compared with the traditional laparotomy, the laparoscopic approach offers the patient the benefit of shorter recovery time.

In our bariatric line of instrumentation, we have available a complete selection of 45 cm long, 5 and 10 mm Ø laparoscopic instruments, with a large selection of grasping, dissecting, biopsy and scissor tips. Rod lens autoclavable 40 cm laparoscopes, extra long trocars and sleeves, the **BRUDER** Anvil grasper for circular staples, 45 cm long fan retractors and the bariatric version of our Endoscopic Fixation Device, which is the **BRUDER** retractor system used for holding laparoscopes, retractors and instruments.

BRUDER instrumentation is known for the highest quality and our fast turn-around repair service for Bruder instruments, as well as for all other brands, is second to none.

***BRUDER** has developed instruments for laparoscopic **gastric bypass and gastric banding**. For these procedures we suggest instruments in two different lengths, 33 cm and 45 cm, to accommodate varying degrees of abdominal wall thickness.*



BRUDER SUGGESTED INSTRUMENTATION FOR BARIATRIC LAPAROSCOPIC SURGERY:

LIST NO. ARO51500	LAPAROSCOPE, AUTOCLAVABLE, 0° VIEW, 10 MM Ø, 40 CM WORKING LENGTH
LIST NO. ARO51501	LAPAROSCOPE, AUTOCLAVABLE, 30° VIEW, 10 MM Ø, 40 CM WORKING LENGTH
LIST NO. ARO10542	RETRACTABLE NEEDLE ELECTRODE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO51510TC	NEEDLE HOLDER, WITH STRAIGHT TIP AND RATCHET AXIAL HANDLE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO51511TC	NEEDLE HOLDER, WITH CURVED TIP AND RATCHET AXIAL HANDLE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO51515	ARTICULATING FAN RETRACTOR, 10 MM Ø, 45 CM WORKING LENGTH

LAPAROSCOPIC INSTRUMENTATION FOR BARIATRIC SURGERY


SUGGESTED INSTRUMENTATION FOR BARIATRIC LAPAROSCOPIC SURGERY (CONTINUED)

LIST NO. ARO51520	ANVIL GRASPER, ROTATABLE, WITH RATCHET HANDLE, FLUSH PORT, 10 MM Ø, 37.5 CM WORKING LENGTH
LIST NO. ARO51147	CLIP APPLYING FORCEPS, MEDIUM-LARGE CLIPS, 10 MM Ø, 37.5 CM WORKING LENGTH
LIST NO. ARO1052715	ASPIRATION NEEDLE, 15 GAUGE TIP, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO10528	SUCTION PROBE WITH CM MARKINGS AND STOPCOCK, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO10533	J-HOOK ELECTRODE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO10535	L-HOOK ELECTRODE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905132B	WAVE GRASPER, ROTATABLE, WITH RATCHET HANDLE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905135B	RETRACTION GRASPER, 2 X 3 ATRAUMATIC TEETH, ROTATABLE, WITH RATCHET HANDLE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905106B	MARYLAND GRASPER/DISSECTOR, ROTATABLE, FULLY INSULATED, WITH RF-CONNECTOR, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905122B	DEBAKEY GRASPER, ROTATABLE, WITH RATCHET HANDLE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905139B	(DORSEY STYLE) FENESTRATED GRASPER, ROTATABLE, WITH RATCHET HANDLE, 5 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905107B	MIXTER GRASPER/DISSECTOR, RIGHT ANGLE, ROTATABLE, FULLY INSULATED WITH RF-CONNECTOR, 10 MM Ø, 45 CM WORKING LENGTH
LIST NO. ARO905165B	METZENBAUM SCISSORS, CURVED BLADE, ROTATABLE, FULLY INSULATED WITH RF-CONNECTOR, 5 MM Ø, 45 CM WORKING LENGTH


For a complete listing of all 33cm and 45cm working length laparoscopic instruments, refer to the "BRUDER Laparoscopic Instruments" data sheet.

The entire line of BRUDER instruments is completely reusable, uncompromising in quality and durability, simple in operation and maintenance and at the most competitive prices.

RECOMMENDED ACCESSORIES FOR BARIATRIC SURGERY SET



LIST NO. ARO100FC	LAPAROSCOPIC TROCAR INCISION WOUND CLOSURE & HERNIA SUTURE RETRIEVER: 3.0 MM DIAMETER, 18 CM LENGTH, SPRING HANDLE DESIGN, POINTED TIP
-----------------------------	--



LIST NO. ARO08801	RETRACT-ROBOT, SINGLE ARM INSTRUMENT: WITH EXTRA-LONG MOUNTING ROD 26", FOR BARIATRIC PATIENTS. (Nathanson Liver Retractor illustrated)
-----------------------------	--

Save and buy direct from the manufacturer.
We design, manufacture & sell the tools the surgeons use.

BRUDER specializes in the design and manufacture of new and innovative instruments, the repair of all competitive products and a fast turn-around time on repairs.



ARROW MEDICAL SUPPLY INC
Telephone Number: (847) 367-7771
Web Site: www.arrowmedical.com

28085 N Ashley Circle Bldg 105

Libertyville IL 600486
Fax Number: (847)367-7751
Email: info@arrowmedical.com